

EXHIBIT 5

Genome Sequencing

Roche Applied Science



- ☐ BIOTECHNICA 2005 Special
- ☐ Genome Sequencer 20 Instrument 3D-flash animation
- ☐ Genome Sequencer 20 System Multimedia Presentation

Literature and References

- ☐ ☐ Sign up for E-Mail Services to stay updated on Genome Sequencer 20 System
- ☐ Download the new brochure
- ☐ Reference: www.nature.com

System description

Revolutionize Whole Genome Sequencing and Assembly

The Genome Sequencer 20 System uses a revolutionary technology, producing more than 20 million raw bases per 5.5-hour run on a single instrument. The software included with the Genome Sequencer 20 System enables mapping or *de novo* assembly for whole genome shotgun sequencing of genomes up to 50 megabases. Many biologically meaningful and complex regions of genomes can be analyzed with this system without the time or cost constraints of current DNA-sequencing methods. The Genome Sequencer 20 System provides an enabling solution for ultra-high-throughput DNA sequencing. In the future, sequencing of complete eukaryotic genomes or high-throughput re-sequencing of human genes within the framework of drug targeting will be possible.

Perform rapid, efficient amplification and sequencing in picoliter format with massive parallelization.

Benefits

- ☐ **Fast**
Sequence more than 20 million bases per 5.5-hour Instrument run.
- ☐ **Cost-Effective**
Benefit from reduced cost per base compared to conventional Sanger technology.
- ☐ **Simple**
Perform sequencing runs with an easy-to-use Instrument requiring minimal steps.
- ☐ **Efficient**
Sequence a typical bacterial genome in days with one person — without cloning and colony picking.
- ☐ **Convenient**
Use the complete system solution — from sample preparation to

data mapping or assembly.




Figure 1: The Genome Sequencer 20 Instrument

DISCLAIMER OF LICENSE: The product is optimized for use in the polymerase chain reaction (PCR) covered by patents owned by F. Hoffmann-La Roche Ltd. No license under these patents to use the PCR Process is conveyed expressly or by implication to the purchaser by the purchase of this product. A license to use the PCR Process for certain research and development activities accompanies the purchase of certain Roche, Applied Biosystems or other licensed suppliers' reagents when used in conjunction with an authorized thermal cycler, or is available from Applied Biosystems. Further information on purchasing licenses may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Dr., Foster City, California 94404, USA.

Restriction of Use. As a condition of sale of this product, purchaser agrees not to use the Instrument to perform less than 10,000 sequencing reactions on a sample or samples in a single sequencing run. Failure to comply with this restriction will result in an infringement of patent rights and other intellectual property rights of seller or third parties and a breach of the terms of the sale of this product.

PICOTITERPLATE is a trademark of 454 Life Sciences, Branford, CT, USA.

Webmaster

© Roche Diagnostics Corporation

Roche Applied Science Home

See our legal Statement